

## 2-(1,2-Epoxy-1,2-dihydrophenyl)acetyl-CoA

## Chemical Properties


CAS No. : 1308875-33-3

Formula: C<sub>29</sub>H<sub>42</sub>N<sub>7</sub>O<sub>18</sub>P<sub>3</sub>S

Molecular Weight: 901.67

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2-(1,2-Epoxy-1,2-dihydrophenyl)acetyl-CoA (2-(1,2-Epoxy-1,2-dihydrophenyl)acetyl-coenzyme A) is a type of acyl coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1091 mL	5.5453 mL	11.0905 mL
5 mM	0.2218 mL	1.1091 mL	2.2181 mL
10 mM	0.1109 mL	0.5545 mL	1.1091 mL
50 mM	0.0222 mL	0.1109 mL	0.2218 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481