


## Plactin B

## Chemical Properties

CAS No. :	174689-42-0
Formula:	C32H52N8O6
Molecular Weight:	644.81
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Plactin B is a cyclic pentapeptide compound that facilitates the dissolution of fibrinolysis proteins.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5508 mL	7.7542 mL	15.5084 mL
5 mM	0.3102 mL	1.5508 mL	3.1017 mL
10 mM	0.1551 mL	0.7754 mL	1.5508 mL
50 mM	0.031 mL	0.1551 mL	0.3102 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481