

PBRM1/SMARCA2,4 binder-1

## Chemical Properties

CAS No. :

Formula: C17H13N3O2

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	PBRM1/SMARCA2,4 binder-1 is an agent that binds to PBRM1, SMARCA2, and SMARCA4, with IC50 values of 0.083 $\mu$ M, 0.153 $\mu$ M, and 0.046 $\mu$ M, respectively.
Targets(IC50)	Epigenetic Reader Domain

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481