

Cy3-Asiatic acid

Chemical Properties

CAS No. :

Formula: C₆₆H₉₇ClN₄O₅

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy3-Asiatic acid is a fluorescent labeling agent formed by the combination of Cy3 and Asiatic acid. It allows for the observation and tracking of Asiatic acid under a fluorescence microscope (Ex/Em = 550/570 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481