


## 3,5,8,11,14-Eicosapentenoyl-CoA

## Chemical Properties

CAS No. :	169337-39-7
Formula:	C41H64N7O17P3S
Molecular Weight:	1051.97
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3,5,8,11,14-Eicosapentenoyl-CoA (3,5,8,11,14-Eicosapentenoyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9506 mL	4.753 mL	9.506 mL
5 mM	0.1901 mL	0.9506 mL	1.9012 mL
10 mM	0.0951 mL	0.4753 mL	0.9506 mL
50 mM	0.019 mL	0.0951 mL	0.1901 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481