


## Ochracenomicin B

## Chemical Properties

CAS No. :	164177-48-4
Formula:	C <sub>19</sub> H <sub>20</sub> O <sub>4</sub>
Molecular Weight:	312.36
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Ochracenomicin B is found in the bacterium <i>Amicolatopsis</i> sp. MJ 950-891, exhibiting no antibacterial activity.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2014 mL	16.0072 mL	32.0143 mL
5 mM	0.6403 mL	3.2014 mL	6.4029 mL
10 mM	0.3201 mL	1.6007 mL	3.2014 mL
50 mM	0.064 mL	0.3201 mL	0.6403 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481