

## 1-Palmitoyl-2-oleoyl-sn-glycerol, 1-Oleoyl-3-palmitoylglycerol (1:1)

### Chemical Properties

CAS No. :

Formula: C74H140O10

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	1-Palmitoyl-2-oleoyl-sn-glycerol,1-Oleoyl-3-palmitoylglycerol (1:1) is a 1:1 mixture of 1-Palmitoyl-2-oleoyl-sn-glycerol and 1-Oleoyl-3-palmitoylglycerol. This compound represents a primary active component of diacylglycerol (DAG) and functions as a hypoxia-inducible factor (HIF) inducer.
Targets(IC50)	HIF/HIF Prolyl-Hydroxylase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481