


(S)-3-Hydroxyoctadecanoyl-CoA

Chemical Properties

CAS No. :	1246036-79-2
Formula:	C39H70N7O18P3S
Molecular Weight:	1049.99
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(S)-3-Hydroxyoctadecanoyl-CoA ((S)-3-Hydroxyoctadecanoyl-coenzyme A) is a type of 3-hydroxyoctadecanoyl-coenzyme A.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9524 mL	4.762 mL	9.5239 mL
5 mM	0.1905 mL	0.9524 mL	1.9048 mL
10 mM	0.0952 mL	0.4762 mL	0.9524 mL
50 mM	0.019 mL	0.0952 mL	0.1905 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481