

SARS-CoV-2-IN-111

Chemical Properties

CAS No. :

Formula: C₃₆H₄₂BrN₃O₃S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	SARS-CoV-2-IN-111 (compound 77) is a highly effective inhibitor of SARS-CoV-2, exhibiting an IC ₅₀ of 0.67 μM. It markedly decreases the presence of SARS-CoV-2 within Vero cells.
Targets(IC ₅₀)	SARS-CoV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481