

PLGA-PEG-PLGA hydrogel (LA:GA 3:1) (20% solution)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PLGA-PEG-PLGA (hydrogel) (LA:GA 3:1) (20% solution) is a biodegradable, thermosensitive hydrogel that undergoes phase transition at specific temperatures. It is widely applicable in drug delivery systems, tissue engineering, wound dressings, and biosensors.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481