

UnyLinker 12 TEA

Chemical Properties

CAS No. :

Formula: C₃₉H₃₅N₁₀O₁₀.1/5C₆H₁₅N

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	UnyLinker 12 TEA is a versatile linker used in the synthesis of oligoribonucleotides.
Targets(IC50)	Nucleoside Antimetabolite/Analog

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481