

## 3(S)-Hydroxy-docosa-7,10,13,16-all-cis-tetraenoyl-CoA

### Chemical Properties

CAS No. :

Formula: C43H70N7O18P3S

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	3(S)-Hydroxy-docosa-7,10,13,16-all-cis-tetraenoyl-CoA (3(S)-Hydroxy-docosa-7,10,13,16-all-cis-tetraenoyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481