


(S)-3-Hydroxyoctanoyl-CoA**Chemical Properties**

CAS No. :	79171-48-5
Formula:	C ₂₉ H ₅₀ N ₇ O ₁₈ P ₃ S
Molecular Weight:	909.73
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.

**Biological Description**

Description	(S)-3-Hydroxyoctanoyl-CoA ((S)-3-Hydroxyoctanoyl-Coenzyme A) is a hydroxy fatty acyl Coenzyme A, formed by the condensation of the thiol group of Coenzyme A with the carboxyl group of (S)-3-hydroxyoctanoic acid.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0992 mL	5.4961 mL	10.9923 mL
5 mM	0.2198 mL	1.0992 mL	2.1985 mL
10 mM	0.1099 mL	0.5496 mL	1.0992 mL
50 mM	0.022 mL	0.1099 mL	0.2198 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481