

E3 Ligase Ligand-linker Conjugate 187

Chemical Properties

CAS No. :

Formula: C₃₀H₄₀N₆O₅S

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	E3LigaseLigand-linker Conjugate 187 is a CRBN E3 ubiquitin ligase ligand-linker conjugate (E3ligaseligand-linker conjugate). This compound is utilized in the synthesis of PROTAC SD-2301.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481