

MP196

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	MP196 is a cationic hexapeptide antibiotic (antibiotic) that targets bacterial cytoplasmic membranes. It exhibits rapid bactericidal action by disrupting membrane integrity, inhibiting cellular respiration, and blocking cell wall synthesis. MP196 holds promise for research on infections caused by drug-resistant bacteria.
Targets(IC50)	Antibacterial,Antibiotic

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481