

Cy7.5 PEG-SH

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy7.5PEG-SH is a fluorescent labeling reagent consisting of CY7.5 dye, polyethylene glycol (PEG), and a thiol group (SH), suitable for molecular conjugation. The SH group is highly reactive and can form stable thioether bonds with various functional groups, such as maleimide (Ex/Em = 770/820 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481