

11-Bromoundecan-1-ol

Chemical Properties

CAS No. : 1611-56-9

Formula: C₁₁H₂₃BrO

Molecular Weight: 251.20

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	11-Bromoundecan-1-ol is a PROTAC linker (PROTAC linker) utilized in the synthesis of PROTACs, such as PROTAC GPX4 degrader-4.
Targets(IC50)	PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9809 mL	19.9045 mL	39.8089 mL
5 mM	0.7962 mL	3.9809 mL	7.9618 mL
10 mM	0.3981 mL	1.9904 mL	3.9809 mL
50 mM	0.0796 mL	0.3981 mL	0.7962 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481