


## GABAA receptor modulator-8

### Chemical Properties

CAS No. : 2889382-34-5  
 Formula: C<sub>20</sub>H<sub>17</sub>FN<sub>4</sub>O  
 Molecular Weight: 348.37  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



### Biological Description

Description	GABAA receptor modulator-8 (10c) is a selective positive modulator of the GABAA receptor that can cross the blood-brain barrier and exhibits anticonvulsant activity.
Targets(IC50)	GABA Receptor
In vivo	The GABAA receptor modulator-8 (10c) demonstrated excellent anticonvulsant activity in the maximal electroshock (MES) model, with an ED50 of 8.51 mg/kg.

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8705 mL	14.3526 mL	28.7051 mL
5 mM	0.5741 mL	2.8705 mL	5.741 mL
10 mM	0.2871 mL	1.4353 mL	2.8705 mL
50 mM	0.0574 mL	0.2871 mL	0.5741 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481