


(3S)-3-Hydroxyadipyl-CoA

Chemical Properties

CAS No. : 1245946-28-4
 Formula: C₂₇H₄₄N₇O₂₀P₃S
 Molecular Weight: 911.66
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	(3S)-3-Hydroxyadipyl-CoA [(3S)-3-Hydroxyadipyl-coenzyme A] is a derivative of coenzyme A.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0969 mL	5.4845 mL	10.969 mL
5 mM	0.2194 mL	1.0969 mL	2.1938 mL
10 mM	0.1097 mL	0.5485 mL	1.0969 mL
50 mM	0.0219 mL	0.1097 mL	0.2194 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481