

L-Serine-S2-PEG-Formic acid

Chemical Properties

CAS No. :

Formula: C₆H₁₁NO₅S₂

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	L-Serine-S2-PEG-Formic acid is a linker used in FGT-4, which acts as a TLR7 agonist and is a chimeric molecule that targets folate receptor beta (FR-β).
Targets(IC50)	Others,PROTAC Linker

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481