


## 2-Methylacetoacetyl-COA

## Chemical Properties

CAS No. :	6712-01-2
Formula:	C <sub>26</sub> H <sub>42</sub> N <sub>7</sub> O <sub>18</sub> P <sub>3</sub> S
Molecular Weight:	865.63
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2-Methylacetoacetyl-CoA (2-Methylacetoacetyl-coenzyme A) is a type of 3-oxoacyl Coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1552 mL	5.7761 mL	11.5523 mL
5 mM	0.231 mL	1.1552 mL	2.3105 mL
10 mM	0.1155 mL	0.5776 mL	1.1552 mL
50 mM	0.0231 mL	0.1155 mL	0.231 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481