

FITC-PEG-Br (MW 3400)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	FITC-PEG-Br (MW 3400) is a fluorescent dye composed of FITC, PEG, and bromine (-Br). The bromine group is an electrophilic active site, suitable for nucleophilic substitution reactions, allowing for subsequent chemical modifications [Ex/Em = 488/525 nm].
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481