

ARCher-142

Chemical Properties

CAS No. :

Formula: C₂₇H₃₁N₆O₆

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ARCher-142 (compound S8) is a potent inhibitor of tankyrase. It selectively binds to ARC4 with an efficacy of 8 μ M.
Targets(IC50)	PARP

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481