


## 3-(m-Hydroxyphenyl)propanoyl-CoA

## Chemical Properties

CAS No. : 86630-22-0  
 Formula: C<sub>30</sub>H<sub>44</sub>N<sub>7</sub>O<sub>18</sub>P<sub>3</sub>S  
 Molecular Weight: 915.69  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-(m-Hydroxyphenyl)propanoyl-CoA (3-(m-Hydroxyphenyl)propanoyl-coenzyme A) is an acyl-coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0921 mL	5.4604 mL	10.9207 mL
5 mM	0.2184 mL	1.0921 mL	2.1841 mL
10 mM	0.1092 mL	0.546 mL	1.0921 mL
50 mM	0.0218 mL	0.1092 mL	0.2184 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481