


## N-Methyl-N-[(3-methyldithio)-1-oxopropyl]-L-alanine

## Chemical Properties

CAS No. :	138148-62-6
Formula:	C <sub>8</sub> H <sub>15</sub> NO <sub>3</sub> S <sub>2</sub>
Molecular Weight:	237.34
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	N-Methyl-N-[(3-methyldithio)-1-oxopropyl]-L-alanine serves as a PEGn linker for antibody-drug conjugates (ADC).
Targets(IC50)	ADC Linker

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2134 mL	21.0668 mL	42.1336 mL
5 mM	0.8427 mL	4.2134 mL	8.4267 mL
10 mM	0.4213 mL	2.1067 mL	4.2134 mL
50 mM	0.0843 mL	0.4213 mL	0.8427 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481