


## Austroside A

### Chemical Properties

CAS No. : 858357-19-4  
 Formula: C<sub>19</sub>H<sub>30</sub>O<sub>9</sub>  
 Molecular Weight: 402.44  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



### Biological Description

Description	Austroside A can be found in <i>Evodia austrosinensis</i> and exhibits antibacterial activity against methicillin-resistant <i>Staphylococcus aureus</i> (MRSA) strains.
Targets(IC50)	Antibacterial

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4848 mL	12.4242 mL	24.8484 mL
5 mM	0.497 mL	2.4848 mL	4.9697 mL
10 mM	0.2485 mL	1.2424 mL	2.4848 mL
50 mM	0.0497 mL	0.2485 mL	0.497 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481