

## 3,5,7-Trioxododecanoyl-CoA

## Chemical Properties


CAS No. : 2484728-87-0

Formula: C33H52N7O20P3S

Molecular Weight: 991.79

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3,5,7-Trioxododecanoyl-CoA (3,5,7-Trioxododecanoyl-coenzyme A) is a type of oxo fatty acyl coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0083 mL	5.0414 mL	10.0828 mL
5 mM	0.2017 mL	1.0083 mL	2.0166 mL
10 mM	0.1008 mL	0.5041 mL	1.0083 mL
50 mM	0.0202 mL	0.1008 mL	0.2017 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481