

Cy5.5 HA MA (MW 50000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

|               |   |
|---------------|---|
| Description   | Cy5.5HA MA (MW 50000) is a methacrylate of hyaluronic acid labeled with CY5.5, widely utilized in the fields of biomarking and imaging, drug delivery, and tissue engineering and regenerative medicine (Ex/Em = 680/710 nm). |
| Targets(IC50) | Others  |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481