

SARS-CoV-2 Mpro ligand 1

Chemical Properties

CAS No. :

Formula: C₂₅H₃₂F₃N₅O₄S

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	SARS-CoV-2 Mpro ligand 1 is a ligand for the target protein PROTAC. It is utilized in the synthesis of [BP-198].
Targets(IC50)	Ligands for Target Protein for PROTAC

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481