


2,4-Dichlorobenzoyl-CoA

Chemical Properties

CAS No. : 150693-18-8
 Formula: C₂₈H₃₈Cl₂N₇O₁₇P₃S
 Molecular Weight: 940.53
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	2,4-Dichlorobenzoyl-CoA (2,4-Dichlorobenzoyl-coenzyme A) is a type of dichlorobenzoyl coenzyme A.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0632 mL	5.3162 mL	10.6323 mL
5 mM	0.2126 mL	1.0632 mL	2.1265 mL
10 mM	0.1063 mL	0.5316 mL	1.0632 mL
50 mM	0.0213 mL	0.1063 mL	0.2126 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481