

LMB-100

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	LMB-100 is a human immunotoxin targeting MSLN/Mesothelin, known for inducing systemic inflammatory response and recruitment of CD8 T cells. It enhances the antitumor effects of anti-PD-1 antibodies in murine models of mesothelioma with human peripheral blood mononuclear cell transplants and syngeneic lung adenocarcinoma expressing mesothelin. LMB-100 is useful for mesothelioma research. Recommended isotype control: Fab Fused with a fragment of Pseudomonas exotoxin A.
Targets(IC50)	Mesothelin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481