

PPAR γ agonist 18

Chemical Properties

CAS No. :

Formula: C₁₆H₁₁BrO₄

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PPAR γ agonist18 is a PPAR γ inhibitor (KD: 3.75 μ M) that suppresses CDK5-mediated phosphorylation of PPAR γ Ser245. It is utilized in research related to insulin resistance.
Targets(IC50)	PPAR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481