

PROTAC PI3K/110 $\beta$  degrader-1

## Chemical Properties

CAS No. : 3070438-85-3


Formula: C<sub>51</sub>H<sub>65</sub>N<sub>9</sub>O<sub>9</sub>S

Molecular Weight: 980.18

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	PROTACPI3K/110 $\beta$ degrader-1 (J-9) acts as a PROTAC degrader specifically targeting PI3K/110 $\beta$ .
Targets(IC50)	PI3K,PROTACs

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0202 mL	5.1011 mL	10.2022 mL
5 mM	0.204 mL	1.0202 mL	2.0404 mL
10 mM	0.102 mL	0.5101 mL	1.0202 mL
50 mM	0.0204 mL	0.102 mL	0.204 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481