


## 3-Hydroxydecanedioyl-CoA

## Chemical Properties

|                   |   |
|-------------------|---|
| CAS No. :         | 138168-88-4   |
| Formula:          | C31H52N7O20P3S  |
| Molecular Weight: | 967.77  |
| Storage:          | Powder: -20°C for 3 years   In solvent: -80°C for 1 year<br>Actual storage temperature shall be subject to the COA. |



## Biological Description

|               |   |
|---------------|---|
| Description   | 3-Hydroxydecanedioyl-CoA (3-Hydroxydecanedioyl-coenzyme A) is a derivative of coenzyme A. |
| Targets(IC50) | Others  |

## Preparing Stock Solutions

|       | 1mg       | 5mg       | 10mg      |
|-------|-----------|-----------|-----------|
| 1 mM  | 1.0333 mL | 5.1665 mL | 10.333 mL |
| 5 mM  | 0.2067 mL | 1.0333 mL | 2.0666 mL |
| 10 mM | 0.1033 mL | 0.5167 mL | 1.0333 mL |
| 50 mM | 0.0207 mL | 0.1033 mL | 0.2067 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481