

Cy7.5 HA MA (MW 10000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from direct sunlight

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy7.5HA MA (MW 10000) is a methacrylate of hyaluronic acid labeled with CY7.5. It is widely utilized in the fields of biomarking and imaging, drug delivery, tissue engineering, and regenerative medicine (Ex/Em = 770/820 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481