

ICG PEG FA (MW 2000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	ICG PEG FA (MW 2000) is a fluorescently labeled biomolecule comprising the Indocyanine green (ICG) dye, polyethylene glycol (PEG), and folic acid (FA). It is widely utilized in drug delivery and biomaterial labeling, with excitation and emission wavelengths of 785/813 nm.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481