


## Lufironil

## Chemical Properties

CAS No. :	128075-79-6
Formula:	C <sub>13</sub> H <sub>19</sub> N <sub>3</sub> O <sub>4</sub>
Molecular Weight:	281.31
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Lufironil (HOE 077) is an inhibitor of prolyl 4-hydroxylase, designed for investigating chronic liver conditions such as alcoholic hepatitis and cirrhosis.
Targets(IC50)	Drug Metabolite

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5548 mL	17.774 mL	35.548 mL
5 mM	0.711 mL	3.5548 mL	7.1096 mL
10 mM	0.3555 mL	1.7774 mL	3.5548 mL
50 mM	0.0711 mL	0.3555 mL	0.711 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481