

## DL-Erythro sphinganine (d20:0)

### Chemical Properties

CAS No. :

Formula: C<sub>40</sub>H<sub>86</sub>N<sub>2</sub>O<sub>4</sub>

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	DL-Erythro sphinganine (d20:0) is a mixture comprising naturally occurring sphinganine (d20:0) and synthetic L-erythro sphinganine (d20:0).
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481