

Cy7 HA MA (MW 100000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy7 HA MA (MW 100000) is a methacrylated hyaluronic acid labeled with CY7, commonly used in bio-labeling and imaging, drug delivery, as well as in tissue engineering and regenerative medicine (Ex/Em = 740/770 nm).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481