


Fluvirucin B3

Chemical Properties

CAS No. :	137120-29-7
Formula:	C ₂₅ H ₄₈ N ₂ O ₅
Molecular Weight:	456.66
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Fluvirucin B3 is a novel antibiotic effective against the influenza A virus.
Targets(IC50)	Antibiotic, Influenza Virus

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1898 mL	10.9491 mL	21.8981 mL
5 mM	0.438 mL	2.1898 mL	4.3796 mL
10 mM	0.219 mL	1.0949 mL	2.1898 mL
50 mM	0.0438 mL	0.219 mL	0.438 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481