

## Astepatide

## Chemical Properties

CAS No. : 2959444-13-2

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Astepatide is an agonist for both glucose-dependent insulintropic polypeptide (GIP) receptors and glucagon-like peptide-1 (GLP-1) receptors, utilized in diabetes research.
Targets(IC50)	Glucagon Receptor,IGF-1R

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481