

(+) -BAY-7081

Chemical Properties

CAS No. :

Formula: C₁₉H₂₇ClN₄O

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(+) -BAY-7081 is an optical isomer of BAY-7081. BAY-7081 is an effective, selective, orally active soluble cyanopyridine-based PDE9A inhibitor with an IC ₅₀ of 15 nM.
Targets(IC ₅₀)	PDE

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481