

ToAP2 TFA

## Chemical Properties

CAS No. :

Formula: C144H235N35O32S.xC2HF3O2

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	ToAP2 TFA is an effective antimicrobial peptide that impacts the membrane permeability of <i>C. albicans</i> and induces changes in the morphology of fungal cells.
Targets(IC50)	Antibiotic,Antifungal

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481