


## 11,12-Didehydroacyl-CoA

## Chemical Properties

CAS No. :	154826-80-9
Formula:	C33H56N7O17P3S
Molecular Weight:	947.82
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	11,12-Didehydroacyl-CoA (11,12-Didehydroacyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0551 mL	5.2753 mL	10.5505 mL
5 mM	0.211 mL	1.0551 mL	2.1101 mL
10 mM	0.1055 mL	0.5275 mL	1.0551 mL
50 mM	0.0211 mL	0.1055 mL	0.211 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481