


2-trans-Eicosenoyl-CoA

Chemical Properties

CAS No. :	71093-69-1
Formula:	C41H72N7O17P3S
Molecular Weight:	1060.03
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-trans-Eicosenoyl-CoA) is a derivative of coenzyme A.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9434 mL	4.7168 mL	9.4337 mL
5 mM	0.1887 mL	0.9434 mL	1.8867 mL
10 mM	0.0943 mL	0.4717 mL	0.9434 mL
50 mM	0.0189 mL	0.0943 mL	0.1887 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481