

[Des-Arg10]-HOE I40 TFA

Chemical Properties

CAS No. :

Formula: C53H77N15O12S.xC2HF3O2

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	[Des-Arg10]-HOE I40 TFA is a potent antagonist of the bradykinin B1 receptor.
Targets(IC50)	Bradykinin Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481