

Macupatide

Chemical Properties

CAS No. : 2923558-68-1

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Description | Macupatide is a gastric inhibitory polypeptide (GIP) receptor agonist, utilized in research related to antidiabetic applications. |
| Targets(IC50) | IGF-1R |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481