


## (E)-2-Benzylidenesuccinyl-CoA

## Chemical Properties

CAS No. : 521064-45-9  
 Formula: C<sub>32</sub>H<sub>44</sub>N<sub>7</sub>O<sub>19</sub>P<sub>3</sub>S  
 Molecular Weight: 955.71  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(E)-2-Benzylidenesuccinyl-CoA ((E)-2-Benzylidenesuccinyl-coenzyme A) is a type of acyl-CoA.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0463 mL	5.2317 mL	10.4634 mL
5 mM	0.2093 mL	1.0463 mL	2.0927 mL
10 mM	0.1046 mL	0.5232 mL	1.0463 mL
50 mM	0.0209 mL	0.1046 mL	0.2093 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481