

Follitropin alfa, Lutropin alfa (2:1)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Follitropin alfa, Lutropin alfa (2:1) is a blend comprising recombinant human follicle-stimulating hormone (Follitropin alfa, FSH) and recombinant human luteinizing hormone (Lutropin alfa, LH) with a 2:1 ratio of FSH to LH. This gonadotropin formulation is used in assisted reproductive technology (ART) to effectively achieve pregnancy while maintaining a good safety profile.
Targets(IC50)	Estrogen Receptor/ERR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481