

E3 Ligase Ligand-linker Conjugate 190

Chemical Properties

CAS No. :

Formula: C29H39N7O4S

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	E3LigaseLigand-linker Conjugate 190 is a conjugate of an E3 ligase and linker, utilized in the synthesis of BP-198.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481